

## SLAMF9 antibody (FITC)

**Cat#: orb13848 (CoA)**

**Cat Number:** orb13848

**Description:** Rabbit Anti-SLAMF9 polyclonal antibody, unconjugated

**Source:** Rabbit

**Clonality:** Polyclonal

**Immunogen:** KLH conjugated synthetic peptide derived from human SLAMF9: 51-150/299

**Purity:** Affinity purification

**Form:** Liquid. Buffer: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**Application:** IF: 1:50-200, Flow-cyt:1:50-200

Excitation Spectrum: 495nm

Emission Spectrum: 519nm

**Storage:** On arrival, store refrigerated at 2°C to 8°C. For long term storage product should be frozen in small aliquots between -20°C to -30°C. Avoid repeated freeze-thaw cycles.

**Note:** Test systems and platform change may cause different results, please contact us if you have any questions.

**Purity Specification:** ≥90% by SDS-PAGE

**Purity Results:** >95% by SDS-PAGE

**Protein/Content Specification:** ≥1.0 mg/ml (specify method)

**Protein/Content Results:** 1 mg/ml by BCA Protein Assay Kit

**Cy3/IgG Ratio Specification:** 3-6

**Cy3/IgG Ratio Results:** ≈4.4